



Seated Support

Instructions



McFee Medical
T E C H N O L O G I E S

800.669.5323

Seated Support

Instructions

The Face Support System, designed and manufactured by McFee Medical Technologies, provided the ultimate in table top or off-the-end of the bed support. The Face Support is adjustable to any height or girth and offers you many daytime activity options. Covered in anti-microbial vinyl, the Face Support features comfortable face padding to support your forehead and cheekbones while allowing you to see.

- 1 Lean your chair forward and lift the back of the seat until the chair is fully open.
- 2 Sit on your chair with your knees on the leg rests and your chest on the chest pad.

Are you comfortable? If so, proceed to Arm Rest Desk Set-Up on the next page. If you aren't comfortable, the seat angle and the height of the chair are both adjustable.

The cable system adjusts the chair height.

- 1 To change the cable setting, pull the chest pad back toward the seat slightly to relax the tension on the cable and remove the "z" clip.
- 2 Insert the "z" clip into a higher or lower hole depending on your comfort. Most patients are comfortable with the "z" clip in the third or fourth hole from the bottom. Higher settings position you higher and more vertically. Lower settings position you lower and more horizontally.
- 3 After you have placed the "z" clip in the appropriate hole, pull back on the seat tube until the cable has no slack. Make sure the "z" clip is parallel to the tube.

The push button on the seat tube adjusts the seat angle and height.

- 1 The back of the seat will be higher than the front with a good basic position.
- 2 Lumbar strain can be reduced by raising the seat height one or two holes.
- 3 To change the seat angle, push in the button on the seat tube while raising or lowering the seat.